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
The gas law apparatus is used in investigations of the well-known special cases of the ideal gas law:

Boyle's law

$$p \cdot V = \text{const.} \quad (\text{for } T, n \text{ constant})$$

Gay-Lussac's law

$$p = \text{const.} \cdot T \quad (\text{for } V, n \text{ constant})$$

Charles' law

$$V = \text{const.} \cdot T \quad (\text{for } p, n \text{ constant})$$

Historically, these gas laws are connected with several physicists. Depending on tradition and language region, each law is known by different names.

Based on the measurements, an experimental value for absolute zero can be found.

Equipment needed

For a complete setup, the gas law apparatus must be combined with the following items, to be purchased separately:

- 007560 Beaker, short, 2 L
- 275010 Immersion heater, 230 V, Schuko plug (or equiv.)
- 064045 Magnetic stirrer (or equiv.)
- 260710 Thermometer (or equiv.)
- 000100 Retort stand base (3 pcs.)
- 000840 Retort stand rod (3 pcs.)
- 002310 Bosshead (3 pcs.)

Demineralized water

- 042300 High vacuum grease

Structure

The apparatus comprises the following parts:

A) A gas volume that can be kept at a well-defined temperature

The volume is divided into a fixed part (aluminium bottle, hosing) and a variable part (plastic syringe). The variable volume is controlled by a spindle with a crank handle. The size is read on the syringe with a resolution of 1 mL.

The temperature of the gas volume is controlled by immersing it into a water bath; it is measured by a standard thermometer.

B) A liquid manometer

The manometer is built around a vertical metre ruler and two glass tubes, connected by a hose. Demineralized water is used as manometer fluid.

The glass tubes can be placed anywhere on the metre stick making it possible to measure pressure differences in a range of about ± 900 mm water column (approx. ± 9000 Pa).

The manometer branch that is open to the atmosphere ends in an overflow container, meaning that this branch can always be placed low while other parameters are being adjusted.

The branch that connects to the gas volume is equipped with a short hose with Luer-Lock, fitting the connection hose.

C) Connection hose with Luer-Lock

The hose is part of a small gas volume that does *not* have the well-defined temperature of the water bath – it must therefore be as small as possible. Therefore an ultra-thin hose is used.

The manometer fluid must be kept out of the hose – otherwise it will be hard to get out again due to the capillary action.

D) Gas amount adjustment system

A short hose goes from the aluminium to a three-way tap. Using an extra syringe, the start pressure can be adjusted somewhat up and down.

The gas volume

As mentioned, the gas we're investigating is divided into a fixed and a variable part:

$$V = V_{\text{FIX}} + V_{\text{VAR}}$$

When analysing the data, V_{FIX} can be considered a known entity (see sections *Calibration of volumes* and *Specifications*) – or it can be treated as a parameter that can be read from a graph.

Calibration of volumes

In section *Specifications*, some measured volumes are given. But these parameters *are not tested* and may vary.

The volume of the bottle can be found by weighing before and after filling it with water with a known temperature. Be careful to fill up the stopper and the pipe stubs as well (use a squeeze bottle for this). All external drops must be dried off.

Hoses etc. can easily be measured with a calliper and a ruler. (With some ingenuity, the water and weighing method may be used here as well.)

Fill the syringe to the 60 mL mark and weigh it. Empty the syringe – but keep the tip filled with water – and weigh it again. This allows the syringe scale to be calibrated.

Note! After using water in hoses, syringe or aluminium bottle, these parts must be dried out *completely* before they can be used for measurements with air! If not, the measurements will be a mixture of air and saturated water vapour.

The manometer

With a thin permanent pen, make a horizontal line at the centre of the glass tube that ends in a Luer-Lock fitting. When the liquid level in this manometer branch is kept at this position in the tube, the fixed volume is kept constant (regardless of where the tube is mounted on the ruler).


Fill the manometer with demineralized water, until the centre of the two glass tubes when these are placed at the same height. Avoid larger air bubbles in the tubes and the hose.

The manometer measures difference pressure rel. to the atmosphere – remember also to measure the barometric pressure (or use the internet).

The pressure difference Δp is given by

$$\Delta p = \rho_W \cdot g \cdot (h_2 - h_1) ,$$

where ρ_W is the density of water, g is the acceleration due to gravity, and h_1 and h_2 are the liquid levels in the two branches.

There may be a small off-set error if the two glass tubes don't have exactly the same inner diameter (due to capillarity). Read h_1 and h_2 with both branches open to the atmosphere and correct the following reading. (Pay attention to the sign.)


The glass tubes of the manometer are fixed to the metre ruler with two plastic clips each. Spare clips are provided.

When the manometer is not used, it should be emptied of the water. Remove the clips from both the tubes and the ruler to preserve their elasticity.

Establishing a possible start pressure

Use the syringe with Luer-Lock and the three-way tap as shown below.

One full syringe volume is more than enough to increase or lower the pressure inside the measurement range of the manometer.

Syringe connected

Syringe disconnected

This changes the amount of substance n of the gas. (We cannot *measure* the actual amount as an independent entity.)

Turn the tap to disconnect the syringe and remove it again.

Airproofing

Hose connections, the gasket under the bottle lid, and the piston of the syringe can get a tiny amount of vacuum grease.

The gas laws

The ideal gas law has the form

$$p \cdot V = n \cdot R \cdot T$$

where p is the pressure of the gas, V is its volume and T is the absolute temperature of the gas.

The amount of substance is n , and R is the gas constant.

In order to establish the ideal gas law, the following special cases has historically played a significant role. The all consider a *fixed* amount of gas.

Boyle's law

$$p \cdot V = \text{const.} \quad (\text{for } T, n \text{ constant})$$

Gay-Lussac's law

$$p = \text{const.} \cdot T \quad (\text{for } V, n \text{ constant})$$

Charles' law

$$V = \text{const.} \cdot T \quad (\text{for } p, n \text{ constant})$$

The three constants are of course not the same.

The gas laws – measurable entities

The gas law apparatus can be used to demonstrate all three gas laws. Typically, you will work with the following rewrites, which can be seen as *linear relationships between appropriately selected (directly measurable) quantities*.

Boyle's law

$$V = \frac{\text{const.}}{p}$$

And when splitting the volume:

$$V_{\text{VAR}} = \text{const.} \cdot \left(\frac{1}{p}\right) - V_{\text{FIX}}$$

(This presents a method to measure V_{FIX} .)

Gay-Lussac's law

$$p = \text{const.} \cdot T$$

Expressed in degree Celsius:

$$p = \text{const.} \cdot t + p_0$$

Charles' law

$$V = \text{const.} \cdot T$$

In degrees Celsius, and splitting the volume:

$$V = \text{const.} \cdot t + V_0$$

$$V_{\text{VAR}} + V_{\text{FIX}} = \text{const.} \cdot t + V_0$$

Significance of the “dead” volume

In most uses of the equipment, the calculations will assume a definite temperature for the gas. But strictly speaking, *a little amount* of gas is not in contact with the water bath. To examine the significance of this small error, we divide the fixed volume into a heated part with volume V_T and temperature T and a much smaller part with volume v and temperature T_A which is not immersed in the water (the “dead” volume).

$$V_{\text{FIX}} = V_T + v$$

Boyle's law

If the experiment is performed at room temperature, all of the gas has the same temperature and the division is unnecessary.

Gay-Lussac's law

From the ideal gas law, the following relation can be derived:

$$p = \left(\frac{\text{const.}}{V_T + v \cdot \frac{T}{T_A}} \right) \cdot T$$

The parenthesis is to a very good approximation a constant, although it will decrease slightly with temperature. Working between 0 °C and 50 °C will only expose deviations from linearity around 1 ‰.

However, if you extrapolate the straight line all the way down to absolute zero, this will typically be determined a few degrees too low.

Charles' law

The gas in the “dead” volume v has constant pressure, temperature and volume – i.e. the amount of gas is a constant as well. This is in other words a completely passive part of the gas. Thus, we can simply subtract this part of the gas volume. (The same conclusion can be reached through calculations based on the ideal gas law.)

Hence the expression changes to:

$$V_{\text{VAR}} + V_T = \text{const.} \cdot t + V_0$$

Experiments

- 132220 Boyle's law
- 132230 Absolute zero (Gay-Lussac's law)
- 132240 Charles' law

Can be downloaded from www.frederiksen.eu

Specifications

Nominal values – not checked!
(Recommended gauged – including the scale of the plastic syringe.)

Heated, fixed volume

Aluminium flask	423 mL
Hose (from flask to syringe)	2,9 mL

Un-heated, fixed volume

Connection hose, thin	1,0 mL
Hose to three-way tap	0,3 mL
Luer-Lock hose + half glass tube	1,6 mL

Variable part

Plastic syringe	0 – 60 mL
-----------------	-----------

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© Frederiksen Scientific A/S