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1975.50 AC

Technical specifications

Weight:	0.45 kg.
Dimension: L x H x D:	160 x 120 x 28 mm.
Temperature range:	0 - 40 °C.
Supply Voltage:	6.3 - 8.0 V via pin 4 and 5
Triggersignal:	Pin 1 and 3 Output: + 3.5V. Input: + 1.0V to + 5.0V.
Rise/Fall time:	< 5 μ s.
Passage time interval:	from 1 ms
Effective width of light beam:	1 mm.

This Photogate Unit is specially designed to be used together with Electronic Counter 2002.50 and the student timer 2002.60. The Photogate Unit is not compatible with Counter type 2001.XX or older.

The Photogate Unit can be used as an optical switch, activated when the light beam is blocked.

The Unit is made of black plastic, with an opening of 90 mm, and a depth of 90 mm. 2 threaded bushings gives the opportunity of mounting the Photogate Unit, using the included mounting rod, in any position wanted.

The Photogate Unit is supplied with a powerful red transmitter diode, and in front of the phototransistor is placed a 1 mm slit. The Photogate Unit is supplied with two LED on the sides of the unit, red for transmitter, and green for receiver, indicating the function.

The Photogate Unit can be connected directly to the DIN sockets of the counter, either one Unit to input A, or two units, one to input A, and one to input B.

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© A/S Søren Frederiksen, Ølgod