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### 134510 Det frie fald – Bestemmelse af $g$

This equipment is designed to be connected to an electronic timer for measuring the time of fall of a free falling steel ball.

The equipment comprises a release mechanism that also acts as a start switch, a target plate also acting as a stop switch, four gold plated steel balls (two  $\varnothing$  12 mm, two  $\varnothing$  16 mm), and a ping-pong ball with steel caps. The ping-pong ball weighs approx. the same as one  $\varnothing$  12 mm steel ball.

A small plastic mirror is also included, It can hang from the rod of the release mechanism to minimise parallax errors in height measurements.


The release mechanism consists of two contact plates that are shorted by the ball used. The ball is held by a spring-loaded neodymium magnet. Good electric contact is obtained by the gold plated surfaces – which should be kept clean and free from grease.

The spring in the release mechanism removes the neodymium magnet fast from the ball. When the ball is released, it drops in a free fall.


When the ball leaves the contact plates, the timer starts. When the ball hits the drop target, electric contact is established there and the timer stops.

The release mechanism is provided with safety sockets for 4 mm plugs, and a  $\varnothing$  10 mm rod for mounting in a stand.


The drop target is provided with safety sockets.


Connection to 200260/200280 Student timer

Release mechanism

Connection to 200250 Electronic Counter

Avoid parallax errors by using the mirror when measuring the height of the fall

### Operations

The setup is shown above.

Prepare the release mechanism by pressing down the button (1) until a click is heard.

Place the steel ball between the contact plates (2).

Use the trigger (3) to release the ball thus starting the timer. When the ball hits the drop target, the timer stops.

### Maintenance and storing

Keep contact plates and balls clean, possibly by polishing with benzine or rubbing alcohol.

For storage, release the spring (press button (3)).

### Experiments

The height of the fall is measured as the distance from *the underside of the ball* (placed in the release mechanism) to *the top surface of the drop target* (depressed).

From the time of fall  $t$  and height of fall  $h$ , the acceleration due to gravity  $g$  can be found, since the following applies:

$$h = \frac{1}{2} \cdot g \cdot t^2$$

A complete lab manual can be downloaded from our website ([www.frederiksen.eu](http://www.frederiksen.eu)):

**134510 The free fall – Measuring  $g$**

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© Frederiksen Scientific A/S