Produktbeskrivelse

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Product descriptions

The switch box is used as an accessory for the airbag bench when acceleration tests are carried out, i.e., that the switch helps to ensure an even start-up of the sled on the air cushion bench during the tests.

The breaker is part of an arrangement where it is inserted between the power supply and the firing mechanism 195210. The breaker box is used to send a start signal to an electronic timer (200270/200285) or a speedgate (197570), at the same time as the current through the coil in the firing mechanism is interrupted. The power supply is a 12 V mains adapter (it is included).

Operation

The switch box is used to send a start signal to an electronic timer (200270/200285) or a speedgate (197570), simultaneously cutting the current through the coil in the firing mechanism.

The power supply is a 12 V mains adapter (it is included).

Connections

In the setup with the air cushion bench, the switch box functions as an electrical switch by connecting it to the electrical switch mechanism and the electrical timer unit. If the setup includes an electric counter in the form of 200270, the setup must follow the principle of Figure 1A.