Betjeningsoversigt

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Operations overview

Preparations: Connect the AC adapter and turn down the *Amplitude* knob.

Switch on by turning the *Waveform* switch from *Off* to the desired waveform (sine, triangle or square). Sine waves will often be the right choice.

Frequency setting: The *Frequency* knob is speed sensitive. Turn quickly to get near the desired frequency and slower to fine tune.

Connection: Connect the *Output* jacks to the speaker, vibrator or other load. The generator has a built-in power amplifier – the output can drive even low impedance loads directly.

Amplitude: Adjust to the appropriate level. (See below how to measure the amplitude.)

Measuring the amplitude

Normal multimeters cannot measure AC voltages across the full frequency range of the function generator. In turn, the output signal amplitude is largely independent of the frequency. Therefore, the amplitude can be measured at a frequency that the multimeter can handle (e.g. 50 Hz), after which the frequency can be set as desired – without touching the *Amplitude* button.

At the highest frequencies, the amplitude will decline slightly - especially for the triangle signal. This poses up to 0.5 dB and can be examined more accurately by using an oscilloscope.

Experiments

The function generator is suitable for most experiments with sound or vibrations and for many measurements on electronic circuits.

Examples of experiments (lab manuals can be found at www.frederiksen.eu):

136210 Rectifier circuits

132810 Resonances in an air column

132860 Standing waves on a string

Specifications

Waveforms:	sine, triangle, square
Distortion (sine):	< 0.25 % (20 Hz - 20 kHz) < 1 % (0.05 Hz - 50 kHz)
Frequency Range:	0.05 Hz to 50.00 kHz
Frequency deviation:	< 0.05 % (for $f > 10$ Hz)
Frequency Drift:	< 0.005 %
Amplitude:	0 to 7.5 V peak (15 V p-p; 5.3 V RMS)
Maximum current:	> 1 A peak (overload protected)
Supply voltage:	12 V (1.5 A) (adapter included)

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© Frederiksen Scientific A/S