
He-Ne Laser 288550

English
manual
at p. 3

2018-04-25 / HS

AC288550

Modulerbar He-Ne laser, klasse 2

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
AC288550

Modulatable He-Ne laser, class 2

This laser emits coherent red light with a wavelength of 632,8 nm.

WARNING!

Class 2 laser (IEC 60825-1)
Do not stare into beam!

The laser beam can burn the retina of the eyes and thus permanently damage the eyesight. Avoid letting direct or even reflected laser beams hit the eyes

The optical power of the laser is up to 1 mW. The light is emitted in a very narrow beam with an intensity high enough to damage the eyesight if the eye is hit.

The laser is provided with a mechanical beam stopper in order to cut off the beam without turning off the laser tube

If a beam with a certain width is needed – e.g. for diffraction experiments – a microscope ocular can be mounted at the front of the laser. This makes a divergent beam that afterwards can be re-focused with a converging lens.

Two sets of steel rods are included for mounting underneath the laser. The longer of the sets fits our optical bench 294600.

The laser is switched on by turning the key switch to a vertical position.

The green LED lights up when power is on.

Careful – high voltage!

If the housing is opened, high voltage circuits are exposed. Internal capacitors will be charged to dangerously high voltages even after the power is switched off.

Modulating

The laser can be amplitude modulated by applying a signal through the BNC socket. This function is primarily for transmission of sound, but frequencies up to 1 MHz can be used if some damping is accepted

Our 489550 Photodetector is suited for use as a receiver (demodulator). It has a built-in loudspeaker.

Typical sources for the audio signal are computers or mobile phones provided with a mini Jack connector. We can provide a cable with a mini Jack plug in one end and two phono plugs in the other (one for each stereo channel). One of these is plugged into a phono-to-BNC adapter.

The BNC connector can be connected directly to e.g. a function generator through a standard BNC cable.

Polarisation

The polarisation plane of the light is not fixed but varies spontaneously. This is especially pronounced during the warm-up phase. If the laser is used with a polarisation sensitive detector the signal will vary as well. This does not imply that the intensity of the light varies.

Additional equipment

- 076886 Objective 10x achromatic
- 288650 Fibre optics Set
- 489550 Photodetector
- 110078 Mini-Jack (stereo) to RCA phono
- 111600 Adapter BNC male / RCA phono female
- 327200 Laser Diffraction Kit

Specifications

Operating voltage: 230 V AC \pm 10 %, 50/60 Hz, 20VA

Fuse: 250 mA, slow blow

Optical power: <1mW

Wavelength: 632.8 nm (red)

Beam divergence (whole angle): 1.8 milliradians

Warm-up period: 20 minutes

Modulation depth: 10-15 %

Input impedance: 1 k Ω

Operating temperature range: 0°C – 40°C

Relative humidity: < 80 %

Dimensions, L x W x H: 460 x 85 x 85 mm

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© Frederiksen Scientific A/S