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15.10.13

3610.65 AC

The power supply no. 3610.65 delivers from 2 to 24 V DC and AC voltage and a current of up to 5 A. The voltage can be selected in increments of 2 V.

The power supply is provided with separate outputs for DC and AC. The two outputs can be used simultaneously. The outputs are protected against overload by a thermal fuse.

The DC output supplies a full-wave rectified voltage without filtering.

Operation

1. Connect the power supply to the mains with the grounded power cord provided.
2. Set the voltage selector to the desired value.
3. Connect the experimental setup to the desired output(s).
4. Turn the unit on. The On/Off button will light up when the device is powered.
5. The voltage can be adjusted while the unit is on.

NB! Please note that the DC and AC outputs are both protected against overload by a thermally activated circuit breaker. If this trips during an experiment, or if there is no voltage on the outputs when the device is turned on: Press the Reset button on the front panel. This resets the circuit breaker. It may of course be necessary to reduce the voltage or reduce the connected load to avoid that the circuit breaker trips again.

NB! The two outlets (DC and AC) must not be connected to each other.

Specifications

Fuse:	4 AT (No. 4090.13)
Power consumption (max.):	180 W
Power out (max.):	120 W
Dimensions: L x W x H:	203 x 225 x 117 mm
Weight:	4.1 kg excl. power cord

DC:

Voltage: 2-24 V (unfiltered), adjustable in 2 V steps
Current: Max. 5 A
Protection: Circuit breaker 6 A

AC:

Voltage: 2-24 volts (RMS), adjustable in 2 V steps
Current: Max. 5 A
Protection: Circuit breaker 6 A

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© A/S Søren Frederiksen, Ølgod