Galvanometer 3815.80

2015.02.23

AC 3815.80

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2015.02.23

AC 3815.80

The galvanometer 3815.80 is a sturdy but sensitive instrument specifically designed for use by students. The instrument is equipped with a zero-centre scale and three current ranges.

Galvanometer applications include zero current measurements used when determining resistance using Wheatstone's bridge. It's easy to balance the bridge by incrementally moving downwards from the 5 mA range to the most sensitive range and here accurately adjust for zero current.

Operation is simple. Connect the circuit between the negative black socket and one of the galvanometer's red sockets, choosing a range higher than the expected current. As a rule of thumb, start with the 5 mA range and move stepwise down to a lower range if necessary.

The current ranges printed immediate to the left of the connectors refer to the full scale reading. For example, a reading of "40" in the 5 mA range indicates that a current of 4.0 mA flows through the instrument.

Frederiksen's student instruments

3815.60 Voltmeter, 3815.70 Ammeter and 3815.80 Galvanometer are a series of basic, rugged analog instruments for student use.

These instruments are built for low voltage use and must not be connected to voltages exceeding 30 V.

Specifications for 3815.80

Precision:	$\pm 5\%$ of full scale
Internal resistance,	$\pm 5\text{ mA}$: ca. $9\ \Omega$
	$\pm 500\ \mu\text{A}$: ca. $90\ \Omega$
	$\pm 50\ \mu\text{A}$: ca. $900\ \Omega$
Dimensions:	158 x 108 x 58 mm
Scale length:	72 mm

The instrument will tolerate up to 1 A overload for at least 10 seconds.

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© A/S Søren Frederiksen, Ølgod