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Student Gas Voltmeter no. 4530.00

Introduction:

For students experiments. The apparatus is made of acrylic which is fairly resistant to chemicals except alcohols and dilutions thereof. Embedded platinum electrodes with terminal bushes for 4 mm plugs. The gas is collected in the two miniature test tubes (incl.).

Electrolysis of a thin dilution of Sulphuric Acid (H_2SO_4) will generate oxygen (O_2) at the positive electrode and hydrogen (H_2) at the negative electrode.

1. Pour a thin dilution of sulphuric acid into the apparatus cup.
The concentration must not exceed 1 M.
2. Fill the two mini test tubes with the diluted acid and place them in the apparatus as illustrated above.
Wear latex gloves.
3. Apply a voltage of max. 5 V DC to the apparatus terminals.

After a while gases will have accumulated in the test tubes. The gases accumulated are hydrogen (H_2) at the negative electrode and oxygen (O_2) at the positive electrode. Hydrogen and oxygen are generated in the ratio 2:1. The result of the electrolysis is a decomposition of the water in its two components hydrogen and oxygen. Likewise electrolysis of other liquids e.g. a thin dilution of hydrochloric acid, salt-water etc. may be done.

Rinse the apparatus thoroughly in running water and dry it carefully before you put it aside.

Trouble shooting:

If the generation of gas is inert – try rinsing the gasvoltmeter in a dilution of ammonia water. Do not use strong dilutions of alcohol.

NB! The gasvoltmeter cannot endure being exposed to ethanol even in weak dilutions.

Technical data:

Cylinder (inner) dia.: 61 mm
Heigth (inner): 56 mm

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© A/S Søren Frederiksen, Ølgod