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


This apparatus demonstrates how an electric current affects a magnetic needle. The large transparent compass rose makes it possible to perform the experiment on an OHP.

The result can be formulated as follows: "Place your right hand with the fingertips pointing in the direction of current, and with the conductor between the magnet and the palm. The north pole of the magnet will then move to the thumb side."

The Danish scientist H. C. Oersted discovered in 1820 this effect after prolonged work with electricity and magnetism. Together with Faraday's discovery of induction in 1831, H. C. Oersted's experiments constituted the initial experimental step towards the complete electromagnetic theory of Maxwell.

Already in 1802 an Italian, G. D. Romagnosi, had observed the deflection of the compass needle, but had not published the results via the usual scientific channels (it was printed in an Italian newspaper).

### Using the apparatus

Turn the apparatus so the conductor is parallel to the magnetic needle. (Because of e.g. metal in the furniture, this is not necessarily the true north-south direction.) Increase the current slowly and observe which way the north pole of the needle is deflected. From 3 to 10 A works convincingly.

As the conductor in this device actually represents a short circuit, it is recommended either to use a power supply *with an adjustable current limiter*, or to *connect an appropriate resistor in series* with the apparatus.

### Recommended equipment

364000 Stabilized Power Supply 0-24 V AC/DC with variable current limiter 0-10 A  
( – or older model: 363000)

– Or a power supply without current limiter in combination with one of these:

420505 Resistor 1 Ω, 5 %, 50 W

420525 Resistor 5 Ω, 5 %, 50 W

**Overskrift**

**Overskrift**

**Specifications**

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© Frederiksen Scientific A/S