Formål med opstillingen

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Purpose of the set-up

With this device, Lenz's law can be demonstrated by letting a magnet roll down an almost vertical 1 m long aluminum rail. Thus, the movement of the magnet can be observed all the way to the end of the rail.

The set consists of:

- 1 x U-profile of aluminum with suspension.
- 2 x cylindrical magnets.
- 2 x non-magnetic "dummy magnets" (cylindrical pieces of stainless steel).

Additional equipment:

- Force sensor.
- Data logger.

Experimental procedure for Lenz's law

1) Set up the U profile so that it is almost vertical. A suitable slope is obtained if there is 5 cm horizontal displacement between the top and bottom. This slope is achieved by allowing the fall ramp to be suspended freely

in the hook as this gives a stable slope of 5 cm. Between the hook and a stand it is possible to place a sensitive force sensor.

2) Let one of the "dummy magnets" roll down the trench - it takes almost half a second to reach the bottom.

3) Then drop one of the right magnets down the trench - this time it takes several seconds before it reaches down.

Theory

If one considers that the sides of the gutter are a closed circuit, then the magnetic flux through the circuit will change when the magnet rolls past. According to Lenz law, a current will be induced in the metal and a magnetic field is formed whose flux through the circuit has opposite sign to the change. In other words, the induced current will counteract the change of the magnetic field.


The magnetic field from the induced current will attempt to repel the magnet as it approaches and retracts it as it exits. The result is that the magnet is braked with a force that is almost proportional to and opposite to the fall rate of magnet. Extension of the experiment

Since the braking force, which affects the magnet (directed upwards), is matched by an equally large and oppositely directed reaction force (downward) on the fall end, the magnet will drop at an even speed (starting speed). (Newton's III law, since the acceleration is 0 m / s^2).

The magnet weighs 4.62 g, which is acceptable as the force which can be measured in the experiment of approx. 0.045 N.

Test run with a power sensor

- 1) Connect a power sensor to the computer or data logger and prepare to take a measurement.
- 2) Connect the U-profile at to the power sensor. Calibrate the sensor and start the data acquisition program.
- 3) Let a "real" magnet roll down the chute until it falls down and stop the data collection.
- 4) Zoom in on the graph (see the example below) over force as a function of time and note that the force is stable as long as the magnet rolls down the U-profile.

(The initial vigorous fluctuations are due to touching the setup at start-up).

Warranty

There are two years of warranty right from the invoice date.

The warranty covers material and production errors. The warranty does not cover equipment that has been mistreated, poorly maintained or incorrectly fitted, and equipment not repaired in our workshop is not covered by the warranty.

Return of defective equipment, which includes Warranty repair at the customer's request, takes place at the customer's expense and risk, and can only be done by agreement with Frederiksen. Unless otherwise agreed with Frederiksen, the freight amount must be prepaid. The equipment must be packed properly. Any damage to the equipment caused by the shipment is not covered by the warranty. Frederiksen pays for returning the equipment after warranty repairs.

© A / S Søren Frederiksen, Oelgod.

This instruction manual may only be copied for internal use at the address for which the appliance is purchased. The guide can also be downloaded from our website, with user login.