Formål med opstillingen

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
Purpose of the set-up

This equipment can be successfully used as a demonstration of how a heat pump works and it is designed specifically for educational use.

The heat pump consists of a compressor, an evaporator and a condenser. The principle behind a heat pump is to transfer heat from the evaporator to the condenser. Here, it is utilized that the phase transition from liquid to gas phase is very energy intensive and the reverse process is similarly energy release.

In order for this process to proceed, an external energy supply provided by the compressor is required.

Content of model

- Compressor
- Evaporator
- Condenser
- Containers for water (not included).
- Thermometer (not included).

Technical data

- Refrigerant: Isobutane (R600a - Freon Free).
- Working pressure, evaporator: 0-3 bar.
- Working pressure, capacitor: 6-15 bar.

Refillable of gas

The system must be refilled with refrigerant at intervals.

1. Unplug the heatpump from power
2. Evacuate the heatpump
3. Fill Isobutan (R600a) until the evaporator manometer shows about 0,9 bar (Approx. 25 g.)
4. Connect the heatpump with power, switch power on
5. Adjust the pressure in the evaporator to 0 bar (working pressure in the evaporator must be 0 bar)
6. After working 5 minutes switch off the heatpump, disconnect the filling station from heatpump und close the valve with the cover.

Theory:

A heat pump raises the temperature of a medium that takes heat from the environment and transfers heat energy for consumption, such as room or water heating. The principle is the same as in a refrigerator. The refrigerator draws the heat out of the closet and delivers it to the room outside. With a heat pump it just the other way around. It takes the heat from the outside and delivers it inside.

Heat pumps use energy to create pressure, and every time they use 1 kWh, 3-4 kWh is created again in the form of usable heat energy.

The heat energy that the heat pump derives from the surroundings is either "free" (eg. temperate air) energy from nature or recycled heat from the home. The "free" energy can, for example, be solar heat stored in the soil, air, groundwater or other heat reservoir (water tank).

The recycling heat can e.g. be excess heat from ventilation systems, wet rooms, roof rooms or stables.

The teaching model provided by Frederiksen has a vaporizer that cools and a condenser that heats. By immersing these in suitably insulated water containers, heat transfer from one heat storage to another can be observed at a higher temperature.

There are differences between how good different types of heat pumps are. The power factor is a measure of the utility or efficiency of the heat pump. The power factor is the ratio of the amount of heat that the heat pump emits (measured in kW) to the power supplied to operate the heat pump (measured in kW).

Functionality


The heat pump works so that heat energy is moved from a cold area to a higher temperature area. These reservoirs are thought to be very large.

Without the use of extra energy, this process would be contrary to Clausius' formulation of the second principle of thermodynamics:

"A process whereby heat energy is transferred from a low temperature area to a higher temperature area is not possible without the application of energy outside the system."

A heat pump is a device that allows energy to be transferred from a cold area to a hot area by supplying energy from outside. Figure 1 shows a principle sketch of this process. Ideally, the applied heat energy ΔQ_V will be equal to the sum of the heat energy ΔQ_K that is moved, plus the energy A supplied from outside:

$$\Delta Q_v = \Delta Q_k + A$$

You can also consider the equation as a power equation, since A is understood as the power output (energy per second). In this ideal case, ignoring energy losses, etc., ΔQ_V will be transferred to an area where is a heat requirement while

providing only (paying for) the energy A . In practice, for example, you will obtain the heat from the outdoor air, from a large water reservoir or similar energy storage with virtually unlimited low temperature heat energy available.

From the laws of thermodynamics one can show that;

$$COP = \Delta Q_k / A = \Delta Q_k / (\Delta Q_v - \Delta Q_c)$$


Where **COP** stands for **Coefficient Of Performance**. In the formula, the counter shows how much heat energy is moved from the cold reservoir. The denominator tells how much energy needs to be sacrificed to move it. The formula thus explains the relationship between the transferred heat energy and the energy for which it is collected (paid for). In heat pumps used, the COP will be lower than the theory here predicts. Theoretical values for the COP of heat pump systems typically range from 2-3 and up to 5-6. Heat loss and other practical conditions mean that realized COP values range from about 1 to 4.

The COP formula can also be used with the temperatures of hot and cold reservoir included (T stated in Kelvin):

$$COP = T_k / (T_v - T_k)$$

Functioning principle of the model

The compressor forces an appropriate fluid from the evaporator to the capacitor, a region of higher pressure and temperature (TV). The fluid is forced further through an expansion valve. The pressure and temperature decrease, and the evaporator is an area at lower temperature (TK). The two heat exchangers shown in the figure allow heat energy from the surroundings to be transferred to and from the heat pump.

An example

A heat pump must transfer heat energy from the outdoor air at a temperature of 10 ° C to a hot water tank at a temperature of 50 ° C => 323 K.

Calculate the maximum theoretical COP.

$$\text{COP} = 283 / (323 - 283) = 283/40 = 7.075$$

Typically, heat pumps used have COP values of about half of this ideal value.

Experiments

Trials with this heat pump can be divided into two types:

1) Quick demonstration attempt:

Demonstrate how a heat pump transfers energy.

1. Place two water-filled containers so that both the evaporator and the capacitor are covered by water. A 1000 ml cup is suitable for this purpose.

2. Turn on the heat pump and observe that the evaporator container is cooled by 5-10

degrees and that the condenser container is heated by 10-20 degrees.

2) Quantitative demonstration experiment:

Make a quantitative experiment over time and with several parameters.

1. The two water tanks should be insulated so that the heat exchange with the environment is avoided at least for a period during the experiment. Use e.g. rock wool supplemented with a flamingo cover.

2. Use temperature sensors to log the temperature evolution over time. (Stir under way in the containers, as there is otherwise a significant temperature gradient in the two heat stores).

3. It is also possible to calculate the COP using EXCEL.

Warranty

There are two years of warranty right from the invoice date.

The warranty covers material and production errors.

The warranty does not cover equipment that has been mistreated, poorly maintained or incorrectly fitted, and equipment not repaired in our workshop is not covered by the warranty.

Return of defective equipment, which includes

Warranty repair at the customer's request, takes place at the customer's expense and risk, and can only be done by agreement with Frederiksen. Unless otherwise agreed with Frederiksen, the freight amount must be prepaid. The equipment must be packed properly.

Any damage to the equipment caused by the shipment is not covered by the warranty. Frederiksen pays for returning the equipment after warranty repairs.

© A / S Søren Frederiksen, Oelgod.

This instruction manual may only be copied for internal use at the address for which the appliance is purchased. The guide can also be downloaded from our website, with user login.