Italiano

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Pre-operation Treatment
of Electrodes

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English

INTRODUCTION

These instructions apply to measuring pH or ORP, reference, combined pH or ORP electrodes.

Each electrode has a bulb protector boot filled with a membrane junction agent which provides an ideal storage for long periods.

REMARK: the electrodes should not be stored for a period longer than 6 months.

PREPARATION

1. Remove the bulb protector boot then immerse the lower end of the electrode into distilled water or a pH buffer for 30 minutes. This operation hydrates the pH bulb and/or the reference junction for optimum performance (not applicable for measuring ORP electrodes).
2. Fill refillable electrodes with the reference solution indicated on the electrode (generally KCl 3M) to a level just below the filling hole(s). Gently shake the electrode downward in the same manner as a clinical thermometer to remove trapped air bubbles.

REMARK: the level must always be kept above the internal element.

pH STANDARDIZATION

Two fresh buffers are required (for high accuracy set the buffer temperature to 25°C).

-pH 7.00 for initial standardization, followed by pH 4.01 or pH 9.22 (choose the nearest to the pH sample) to span the electrode (do not forget to rinse the electrode(s) between the two measurements).

ORP STANDARDIZATION

One or two ORP standards are recommended to check that the electrode is working properly.

ELECTRODE CLEANING

Electrodes are susceptible to coating by many substances and the response time can deteriorate dramatically. If they are mechanically intact they can often be restored to normal performance by one of the following procedures.

-**salt deposits:** HCl 0.1M (5 min) then NaOH 0.1M (5 min)

-**Oil/grease films:** detergent and water

-**Clogged reference junction:** KCl 0.1M (60-80°C during 10 min)

-**coated/abraded (for platinum ORP):** aqua regia¹ HCl concentrated/HNO₃ concentrated (1/1).

SPECIFICATIONS

They are indicated on each electrode.

¹ This procedure uses hazardous chemicals and should only be performed by a qualified person.

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