Plasts egenskaber (bestrålet), 582010

Produktbeskrivelse

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Plastic's capabilities (Radiated), 582010

Product descriptions

The set can be advantageously used for studies of the properties of plastics under various influences. Inspiration for this can be found in the book "Chemical production - and your everyday life" published by the publisher Malling Beck.

By using the kit, the following different properties of plastic can be examined:

- Electron irradiation of plastic gives higher flexibility.
- Certain types of plastic contract when exposed to heat.
- Hydrophilic property of certain types of plastic.
- Property for PVC with plasticizer.

Operation

The following experiments can be supported by the kit:

Dried sponge is water absorbent.

- Shrink wrap shrinks.
- Irradiated and non-irradiated PLASTIC
- Plasticizer in PVC

Experiment 1 - Dry sponge is water absorbent.

Necessary equipment:

- 3 pieces of washing sponge.

- 3 beakers.
- Rubbing alcohol.
- Water.
- Cooking oil.

Procedure

- Pour the liquids into separate beakers and put a small piece of dried sponge in each glass. Watch what happens to the three pieces of sponge.
- As the washing sponges are hydrophilic ("water-loving") will
- The sponge in water absorbs this and grows. Compare with the other two mushrooms (nothing will happen to them).

Test 2: Shrink plastic shrinks.

Necessary equipment:

- A piece of shrink wrap (approx. 5 cm).
- Wooden pencil, or wooden stick of the same diameter.
- Candles.

Course of action

- Place the shrink-flex hose halfway on the pencil/
- the wooden stick and hold it carefully over the candle flame so that it heats up slowly.
- Turn the pencil/wood stick slowly so that the heat is evenly distributed over the entire shrink wrap.
- Observe what happens to the shrink wrap around the pencil and what happens to the shrink wrap outside the pencil/stick.

Additional experiment: Try pouring boiling water into a plastic soda bottle and observe what happens.

Experiment 3, Irradiated and non-irradiated plastic - heating of PE balls

Theory:

When PE is irradiated, cross-links are formed between the individual polymer chains. This changes the plastic from having an almost linear structure to having a net

structure. This causes the plastic to become more cohesive when heated. The unirradiated PE ball will burst when pulled, while the irradiated PE rod will stick together. Similarly, the unirradiated PE ball will burst and flatten after heating, whereas the irradiated PE ball will retain its shape.

Necessary equipment:

- PE balls, irradiated and non-irradiated.
- Beaker, 500 ml.
- Cooking oil.
- Tripod.
- Ceramic wire mesh.
- Thermometer, -20 – 360 °C.
- Bunsen burner.

Course of action:

- Pour approx. 250 ml of cooking oil in the beaker and heat it to around 150 °C.
- Put the two PE balls into the cooking oil and watch what happens.
- Carefully remove the hot balls from the hot oil, place them on a piece of paper towel and let them cool.
- Describe how the two balls look and feel after cooling

Experiment 4 - Plasticizer in PVC**Theory:**

PVC is an excellent plastic material for many things, and is in itself rigid, which is why plasticisers, the so-called phthalates, are added (under suspicion of disturbing the hormonal balance in living organisms).

Different solvents affect PVC and the added plasticizer in different ways, so you must carefully consider what it will be used for before using PVC - as with all plastic materials.

Necessary equipment:

- PVC hose.
- 4 glass containers with lids.
- Water.
- Acetone.
- Methanol.

- Propanol.

Course of action:

- Cut four small pieces (5-10 cm) of plastic tubing and place them in separate containers with each of the above liquids.
- Leave them for a few days and then pick them up.
- Observe what has happened to the plastic.

The experiment may possibly be supplemented with other liquids.