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```

on start
  radio set group 1

on button A pressed
  motor 1 on direction forward speed 30

on button B pressed
  turn off motor 1
  
```

```

on radio received receivedString
  if receivedString = 'Dry' then
    show icon
    motor 1 on direction forward speed 20
    pause (ms) 100
    play sound
      start frequency 5000
      end frequency 0
      duration 500
      until done
    pause (ms) 350
    turn off motor 1
    clear screen
  
```

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The 3 V water pump is perfect for creating your own self watering greenhouse, or water fountain project.

Necessary accessories:

Simple Water Fountain:

- 663039 water pump
- 038530 silicone tube
- on/off switch
- 641320 wire connector
- extra wire

Create your own simple water fountain by connecting the water pump directly to a 3 to 4.5 V power source. We recommend a 4.5 V battery or the universal power supply (355035) Use a wire connector to extend the wires of the water pump (soldering the small wires can be difficult), so that the power supply can be connected away from the water. Add a switch into the circuit. The pump will operate at full capacity as long as the switch is turned on. If you are interested in regulating the speed of the pump, see the instructions for the self watering greenhouse below.

Self Watering Greenhouse


- 663039 water pump
- 038530 silicone tube
- 663013 motor controller
- 641320 wire connector

- 663002/663020 microbit
- 663024 soil moisture probe
- extra wire

In your self watering greenhouse, you will be able to determine the moisture of the greenhouse/plant using a soil moisture probe. You will use the micro:bit to read the value of the probe, in order to determine if you need to turn on the water pump. The motor controller is necessary as its outputs are designed to deliver sufficient current to drive the water pump. You might be tempted to connect the water pump directly to P0, P1, or P2 on the microbit, as they can deliver a 3 V signal, but they cannot deliver more than 15 mA. That is why we will be using a motor controller.


Connect the red motor wire to P8 on the motor controller and the black motor wire to P12 on the motor controller. To extend the black and red motor wires, connect extension wires using the 641320 wire connector.

In the code example below, we used two micro:bits. To make the project more interesting, we send a message from the soil moisture system to the water pump system when it is time to water.

Soil moisture sensor code:

Water pump code, including manual start and stop of pump. The speed shown below is 30% of maximum pumping speed.


```

on start
  radio set group 1

forever
  set Moisture to analog read pin P1
  show number Moisture
  pause (ms) 500
  if Moisture <= 500 then
    radio send string Dry
  
```

Water pump code, continued. Receiver code to water when soil sensor is dry.


```

on start
  radio set group 1

on button A pressed
  motor 1 on direction forward speed 30

on button B pressed
  turn off motor 1
  
```

```

on radio received receivedString
  if receivedString = Dry then
    show icon
    motor 1 on direction forward speed 20
    pause (ms) 100
    play sound
      start frequency 5000
      end frequency 0
      duration 500
      until done
    pause (ms) 350
    turn off motor 1
    clear screen
  
```

Alternative project:

If your watering greenhouse does not need to regulate the speed of the water pump, you can replace the motor controller with a relay or switch.

Another important feature of a motor controller is that it can give you negative voltage, for example to make a motor go in reverse, which we will not be doing here.