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Frederiksen Scientific ready-to-use kit for qualitative production of biogas.

The kit includes 3 bottles with lids, acrylic tubes, silicone hoses and hose clamps, as well as a product guide and an experiment manual for a simple laboratory experiment with biogas production from decomposing organic material.

The gas production can be tested by flame test or collected and tested on a gas chromatograph.

Application

The biogas kit can be used in three different ways.

1. Testing if biogas can be produced from a selected organic material. Here, the focus is on gas production, rather than the amount of gas produced. Gas production is tested by flame test.
2. Comparing gas production from different organic materials and establishing which material produces the largest amount of gas during the given test time. The gas quantity is measured on **Bottle 3**.
3. Testing gas production from the same organic material under variable conditions. Varying parameters could be: temperature, the amount of starter material – if such is used, the surface area of the organic material, stirring or not stirring, etc. The gas quantity is measured on **Bottle 3**.

Materials

076630 Frederiksen Scientific Biogas kit

085807 Lab coats

086047 Disposable gloves

085030 Safety glasses

300 g organic material - finely shredded

Possibly starter material

Water

064066 / 065210 Magnetic stirrer preferably in combination with heat source

065520 Magnets

Various lab equipment such as: beakers, funnels and other glass containers.

Organic matter and starter material

Biogas can be produced from many different types of organic material – it could be household waste or waste from food industry and agriculture. Most organic materials already have bacteria capable of decomposition on their surface. Some of these bacteria will, under anaerobic conditions, produce methane (CH_4) during the decomposition process.

To ensure a quick start-up and methane production, it may be an idea to get or buy a small amount of starter material from a local biogas plant. If you don't have access to it, pigs faeces can be used as starter material.

In terms of quantity, you don't need more than approx. 25 ml of starter material for each biogas kit.

If you want to use specific organic materials for gas production, it may be necessary to acquire specialized bacterial cultures.

Preparation:

Start out by finely chopping the organic matter from which you want to extract biogas. A chopping ensures a larger surface area from which the decomposition can take place.

If possible heat the experiment to 35-37 °C. Higher temperatures will speed up the biochemical reactions by increasing the metabolism and division rate of the methane-producing bacteria.

Setup

Start by filling the 3 bottles as described below.

Bottle 1: Add 300 g of finely chopped organic material and add, if possibly a little starter material – approx. 25 ml. Then fill the bottle to above the 400 ml mark with water. Put a magnet in the bottle and put it on a magnetic stirrer, stirring slowly. If the magnetic stirrer has a heating option set temperature to 35-37 °C.

Biogas is produced in **Bottle 1**.

Bottle 2: Fill the bottle with water. The less atmospheric air above the water, the more concentrated the gas that collects in the bottle will be.

The biogas is stored in **Bottle 2**.

Bottle 3: Kept empty from the start. This bottle acts as a 'pressure valve' and a measuring device.

It is in **Bottle 3** the liquid the biogas displaces from **Bottle 2** ends up. The more liquid ending up in **Bottle 3**, the more biogas is produced in **Bottle 2**. The liquid in **Bottle 3** thus becomes an estimation of the total biogas production in **Bottle 2**, and makes it possible to compare gas production from different types of organic material, or gas production from the same type of organic material tested under varying conditions.

The prepared bottles can be seen above.

When the three bottles are ready, assemble the 'biogas plant' as shown in C below. Corks, tubes and pipes must fit snugly together to avoid the produced gas leaking out during the experiment.

Note that the hose clamps are threaded on the three connecting tubes, but **NOT** closed. The hose clamps are

ONLY closed when the produced gas from **Bottle 2** is to be tested or moved out of the bottle.

If you have a heating option available you can run the experiment at a temperature of approx. 35-37 °C. Otherwise, the test is carried out in a fume cupboard in a warm room.

The experiment lasts for 7 – 21 days. When the water level in **Bottle 3** has not changed for a few days, gas production has ceased and the experiment can be stopped.

Test results

If the water level in **Bottle 3** rises and the water level in **Bottle 2** falls, gas is produced in **Bottle 1**.

The gas produced in the experiment is generally a mixture of two gases: carbon dioxide (CO₂) and methane (CH₄). In order for methane to form, the environment in **Bottle 1** must be anaerobic - i.e. free of oxygen. When the experiment is started, O₂ will be present in **Bottle 1** and CO₂ will hence be formed via bacterial respiration. When the oxygen in **Bottle 1** is used up, methane production begins.

If the gas that has been produced contains a large CO₂ fraction, it cannot burn. Therefore, you can test whether sufficient methane has been produced in **Bottle 1** by closing the two hose clamps as close to **Bottle 2** as possible, disconnecting **Bottles 1** and **3** and then testing the gas in **Bottle 2**.

Remember to write down the water level in **Bottle 3**.

The easiest way to test the gas is a qualitative test where you light a match, remove the silicone tube previously attaching **Bottle 2** to **Bottle 1** and quickly hold the match towards the mouth of the acrylic tube. If the flame is extinguished the CO₂ content is too high. If the gas burns with a blue flame, methane and thus biogas has been produced.

If you want to collect the gas for other tests, possibly gas chromatography, the gas can be transferred from **Bottle 2** into a urine bag via a disposable syringe (100 ml).

Detach **Bottle 2** from the other bottles as described above. Fill the syringe with water, press the air out of the hose that originally connected **Bottle 2** with **Bottle 1** – and connect it to the syringe. Attach the urine bag to the silicone tube that connected **Bottle 1** to **Bottle 2**. Again squeeze the air out of the tube before you connect the urine bag, to minimize contamination by atmospheric air.

Now press the water from the syringe into **Bottle 2** – you will need more than the 100 ml of water. The water filling up **Bottle 2** will transfer the gas from **Bottle 2** to the urine bag.

Security measures

- Use gloves when handling the organic material – especially if it is not fresh.
- Use gloves when handling the starter material or if pig faeces is used.
- Keep an eye on the water level in **Bottle 3** and make sure to drain **Bottle 2** of gas when the water level in **Bottle 3** becomes high and the water level in **Bottle 2** becomes low.
- Take care when flame testing the gas. Make sure to wear safety glasses and keep the tube opening away from people.
- Be aware that the biogas is flammable, so if it is collected in large quantities, it must be handled with care and kept away from open flames.
- Use gloves when the experiment is ended and especially when **Bottle 1** is washed.
- The contents of **Bottle 1** must be incinerated, or buried in a field or the likes.