

2021-11-03 / MMK

Productmanual AE\_077025/077026 p. 1/1



# MANUAL FOR STEREO MICROSCOPES: Model 30B: 10x and 30x magnification Model 30C: 20x and 40x magnification

- If the black/white stage plate is mounted then please replace it with the glass stage plate. The frosted side of the glass stage plate should face down.
- Remove the cover from the underside of the objective tube (not always mounted).
- Put the mains cord in a socket and turn on the stereo microscope. The model has two illumination sources; top light (incident light) and bottom light (transmitted light).

### General use of the stereo microscope:

- Place an object under the objective revolving tube, and use the focusing knob to adjust the sharpness.
- If a sharp image cannot be attained due to the height of the object, then carefully release the lock screw on the vertical post and find the correct height. Remember to lock the lock screw again. The support collar also needs to be locked.
- The magnification can be changed by turning the vertical revolving objective tube a half turn until a click is observed.
- Place the black/white stage plate on top of the glass stage plate, if you want to see an object on a black or white background. Incident light is used for this purpose.

- The glass stage plate should not be removed. Using the black/ white plate on top of the glass plate prevents the plastic plate from damage if the transmitted light is not turned off.
- Interpupillary distance is adjusted on the stereo microscope by moving the eyepiece tubes/prism housings. Look into the eyepieces and move the prism housings (light colour) closer together or further apart. The interpupillary distance is adjusted when a single circular field is observed when viewing with both eyes.
- The stereo microscope is then adjusted to your sight to achieve a good and sharp image. Look through the right eyepiece with your right eye. Use the focusing knob to get a sharp image of the specimen. Then look with your left eye through the left eyepiece and adjust the sharpness with the dioptre ring on the eyepiece tube. The adjustment is finished, when the image seen with the left eye is sharp.

## Adjustment of the height:

- The stereomicroscope is mounted on a pole stand. As mentioned above, this makes it possible to see large (and high) objects under it. Please be aware that the stereo microscope gets unstable, if the optical head is turned away from the main stand. It is important to fasten the stereomicroscope to a table if it is used in such a manor.
- The friction in the focussing adjustment can be altered with the enclosed special key if the stereomicroscope cannot hold the focus because it drops down. It is also used if the focusing knob is very tight to turn. The tip of the key is placed in one of the four holes in the axle inside one of the focusing knobs. It is then possible to tighten or loosen the ring that is placed around the axle, and with it adjusting the friction.

### Bulb:

Incident light (top): Halogen 12V/10W Transmitted light (bottom): Pinol Bulb 12V/10W

## Accessories (not included):

#### Eyepiece 20x.

Remember to turn the power switch off and to disconnect the power cord from the main supply, when the bulb is changed.

The longevity of the stereomicroscope can be improved if it is protected against dust, humidity, high temperatures and shocks.