# Micro Capture Software User Manual Ver7.2

# CONTENT

CHAPTER 1: SOFTWARE INSTALLATION	- 1 -
1.1 Environmental Requirements for Software Installation	1 -
1.2. Installation	1 -
CHAPTER 2: SOFTWARE INSTRUCTIONS	6 -
2.1 Software Operation Interface and Instruction	6 -
2.2 Manual Panel	7 -
CHAPTER 3: FUNCTION INTRODUCTION	- 8 -
3.1 Shortcut toolbar	- 8 -
3.2 Working mode	11 -
3.3 Camera control	- 13 -
3.4 Measuring tool	16 -
3.5 Drawing Tools	- 18 -
3.6 Annotation Tools	19 -
3.7 Eraser	- 20 -
3.8 Image Filters	21 -
3.9 Calibration and Measurement	- 22 -
3.10 Thumbnail display area	27 -

# **Chapter 1: Software Installation**

#### **1.1 Environmental Requirements for Software Installation**

Compatible system: XP, Win7, Win8, Win8.1, Win10

Note: Partial function of this software relies on the installation of the following software: Microsoft Word 2003 or later version (for exporting report to WORD) Microsoft Excel 2003 or later version (for exporting report to EXCEL) Microsoft Outlook 2003 or later version (for exporting images to email) Internet Explorer7.0 or version above (for wireless connection to PC)

### **1.2. Installation**

Step 1.Activate the computer and insert the installation disk into the DVD.


Step2.Double click the file"Micro Capture Install" under the disk catalogue. When the following window pops up, please select the installation guide language and press the button"

Select Set	up Language
12	Select the language to use during the installation:
	English
	OK Cancel


Note: The installation guide language is independent of the software interface language. If you need to change the software interface language, please refer to 2.2 - Modification of Language Function.

Step 3. When the window" Welcome to use Micro Capture Installation Guide" pops up, please

click the button "



Step 4. When the window "Please select the target location" pops up. click "Browse" to choose the file path for installation (usually default), and then click " Next > " to continue.



Step 5. When following window appears, select the shortcut you need, and then click the button



追 Setup - Micro Capture	
Select Additional Tasks Which additional tasks should be performed?	
Select the additional tasks you would like Setup to perform while installin Capture, then click Next.	g Micro
Additional icons:	
Create a desktop icon	
Create a Quick Launch icon	
< Back Next >	Cancel

Step 6.Confirm information like "the target location" of software installation. After


confirmation, please click the button

Setup - Micro Capture	
Ready to Install Setup is now ready to begin installing Micro Capture on you	ır computer.
Click Install to continue with the installation, or click Back if change any settings.	you want to review or
Destination location: C:\Program Files (x86)\Micro Capture Start Menu folder: Micro Capture Additional tasks: Additional icons: Create a desktop icon	*
	*
< <u>B</u> adk	Install

Step 7. An installation progress bar will appear during the software installation procedure. The software installation necessitates the registration of some runtime libraries, which might be intercepted by some antivirus software. If such interception pops up, please allow the procedure to continue.

Step 8. An installation progress bar will appear during the software installation procedure. When the following window appears, the software has been installed successfully. Then click the button

**Einish** to complete the installation. After closing the window, check the icon **P** on your desk. If you fail to find it, please find it from the path "Start" $\rightarrow$ "All programs".




Note: This software supports connecting the camera to PC through wired or wireless way.

Wired connection: just connect the camera to USB port of PC. If it is a WiFi camera, please dial the hardware switch of WiFi camera to USB working state in advance.

#### Wireless connection: only valid for WiFi camera.

- 1. Please switch the hardware switch of WiFi camera to WiFi working state in advance
- 2. Connect your computer to the WiFi camera, please confirm that your computer has a wireless network in advance. Since the default of our camera is 5G WiFi IEE 802.11ac protocol, please confirm in advance that your computer wireless network supports this protocol. If not, please refer to the hardware specification of the camera and switch the camera to 2.4G WiFi IEE 802.11b before connecting with the PC.
- 3. After starting the WiFi camera, search for WiFi name (WiFi camera xxxx) in the network connection. The default connection password is 12345678




4. Run the software

Select the "<sup>O</sup>" button in the shortcut toolbar, and then select" WiFiCamera "from the"

WiFi Camera "selection in the upper right corner of the software. Wait a moment for the PC to connect to the WiFi camera.

# **Chapter 2: Software Instructions**

# **2.1 Software Operation Interface and Instruction**



Three function bars Video/document preview window

Thumbnail display window

- Working mode
- Camera control
- Measuring tool

# 2.2 Manual Panel

File	Edit	View	Acquire	Process	Measure	Report	Window	Setting
------	------	------	---------	---------	---------	--------	--------	---------


#### The menu Functions are as follows:

Open Image     Save Image     Save Image As     Import Parameter File     Export Parameter File     Exit	▷       Select         □       Drag         ↓       Undo         →       Redo         □       Bring Forward         □       Send Backward         □       Copy         ↓       Crop         □       Paste         □       Delete	<ul> <li>Zoom In</li> <li>Zoom Out</li> <li>Zoom 100%</li> <li>Fit to Screen</li> <li>Full Screen</li> <li>Counterclockwise</li> <li>Clockwise</li> <li>Flip</li> <li>Mirror</li> </ul>	<ul> <li>Preview</li> <li>Pause Preview</li> <li>Capture</li> <li>Timelapse Capture</li> <li>Record</li> <li>Setup Device</li> </ul>
----------------------------------------------------------------------------------------------------------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	--------------------------------------------------------------------------------------------------------------------------------------

Proce	ess Measure Report W	Mea	sure Report Window Setti	0	rt Window Setting	~	Common Tools
$\nabla$	Image Filters	4	The second se	Керо	rt Window Setting	~	View Tools
p	image riners	3/c	Calibration	r ha		~	Editing Tools
G	Drawing Tools	th	Reticle Crosshair		Export to Word	~	Screen Splitting Tools
FG.	Annotation Tools	4	Neucle crossilaii	2000	315 225 2	~	Image Process Tools
-	Francis	22	Measurement Tools	R	Export to Excel	~	Left Toolbar
8	Eraser	<b>P</b> /V	and the second	5		Y	Right Toolbar

Setting	1	
Language 🛛 🕨		中文简体
Picture save path	<b>~</b>	English

# **Chapter 3: Function Introduction**



# 3.1 Shortcut toolbar

	۵۵ 🔍 ایق 🕑 🔂 🕕 🔄 🗈	s - ,0	폐X™ (C ◀ ▲ ད ▻ » ⊕ ⊑ ▣ ☞ ╕ ◻ ⊄ □ □ ⊞ ∀ ☞ ৫ ⊘
D	Open Image	:	Open the image from disk.
B	Save the image	:	Save the image to disk.
	Export to Execl	:	Export the images to the Excel
	Export to Word	:	Export the image to Word
Û	Delete the image	:	Deletes the current image from disk.
igodoldoldoldoldoldoldoldoldoldoldoldoldol	Preview	:	Live-view from camera

	Pause	:	Pause live-view from camera
67	Take picture	:	Capture image to disk
$\odot$	Time-lapse snapshot	:	Time-lapse snapshot to disk
	Record	:	Record to disk
Ð	Zoom in	:	Increase magnification
Q	Zoom out	:	Decrease magnification
(11)	Actual size	:	Display the current image at its actual size
X	Adapt to screen size	:	Display the current image adapt to to screen size
2	Full screen	:	Display the current image at full screen
C	Counterclockwise	:	Rotate the current image 90 degrees counterclockwise.
	rotation		
С	Clockwise rotation	:	Rotate the current image 90 degrees clockwise
<b>V</b>	Flip	:	Flip the current image
Δ	Mirror	:	Mirror the current image
←	Undo	:	Undo the last operation, valid only for the edit operation.
$\rightarrow$	Redo	:	Redo the last action, valid only for editing
$\triangleright$	Select the object	:	Select an object on the screen, only valid for editing objects on
			the screen
Ð	Drag	:	Drag video screen in preview area
9	Flip the object up to	:	Place the selected object on a higher layer
	one more layer		
<b>G</b>	Flip the object down to	:	Place the selected object in the next layer
	one more layer		
∎়	Delete the object	:	Delete the selected object

Ţ	Duplicate object	:	Copy the selected object
Ū	Paste object	:	Paste an object that has been cut or copied.
Ц.	Cut object	:	Cut the selected object
	One screen	:	One screen mode
	Two screens	:	Two screens mode
⊞	Four screens	:	Four screens mode
$\mathbf{Y}$	Special effects	:	Special effects tools
	Drawing tools	:	Drawing tools
ľ	Annotation tools	:	Annotation tools
$\diamond$	Erase tools	:	Erase tools
≣	Detailed information	:	The gallery is arranged in detail
~	Thumbnail	:	The gallery is arranged in thumbnail
<b>9</b>	Slide	:	Display slides
٥, U	ISB Camera	▼	]

This item displays devices currently in use. When the software searches for multiple video devices, multiple device names can be displayed, and devices can be selected, as shown in the figure



# 3.2 Working mode

Microscope Types
Biological Microscope
Industrial Microscope
🗾 Light Source
LED
Low-brightness Halogen
Lamp O High-brightness
Halogen Lamp
Preview Resolution
2592×1944 💌
Capture Resolution
2592×1944 🛛
Capture Format
.jpg 🛛
Microscope Types
Biological Microscope
Industrial Microscope

Microscope options:

Please select your microscope type for better image display.

Light Source
LED
<ul> <li>Low-brightness Halogen Lamp</li> </ul>
<ul> <li>High-brightness</li> <li>Halogen Lamp</li> </ul>

Light source options:

According to the actual situation, select the corresponding light source to get better image display effect.

 $\star$  The above functions are only available for specific models. If there is no display, it means the product does not have this function.

□       Preview Resolution         2592×1944        ▼         1020×1080       -         2592×1944       -         1020×1080       -         2592×1944       -         1024×768       -         3640×480       -	Preview options	:	You can select any of the preview resolutions in this drop-down box to preview as you wish.Note: the higher the preview resolution, the lower the frame rate.When you preview at 1280×720 or 1920×1080, save only 1280×720 and 1920×1080.
Capture Resolution 2592×1944  ▼ 2592×1944 1024×768 640×480	Save resolution	:	This option is used to select the resolution of the image taken.The preview resolution is in proportion to the saved resolution.When preview resolution ratio is 4:3, save resolution only supports 4:3 resolution selection.The default maximum saving resolution is 2592×1944;When preview resolution is 16:9, save resolution only supports 16:9 resolution scanning, default 16:9 save resolution is 1280×720.
Capture Format .jpg V .jpg .bmp .png .tiff .pdf	Save image format	:	There are four save formats: PDF, BMP, JPG and PNG. The default format is JPG.Depending on your requirements, you can choose any of these save formats.
Locked: 🔎 Unlocked: 🛅	Locked and unlocked	:	When it is locked, the preview resolution is consistent with the photo resolution, and the preview screen displays the same as the photo taken. When in the unlocked state, the preview resolution is inconsistent with the photo resolution, and the preview screen may be inconsistent with the photo taken.

### 3.3 Camera control

According to your needs, you can adjust the corresponding parameters to achieve the desired preview and image effect.

As shown in the figure:

Ê	-0-0-	×
-¤́- Brigl	ntness	
Ť	+	45
🖄 🔽 Ex		138
🛞 Area		
LAE		
🖒 Singl	le Exposure	
EV		
🖓 🔽 🗤	/hite Balance	
R	-0	126
В	-0+	
Area	White Balan	
Ewe		
	le White Bala	nce
WBS		
• Cont		
~	+	40
ီး Satur	ration +	64
7	· ·	
🚳 Shar		
(	· ·	50
Powe		
50Hz	∞  ▼	
	eve Settings	



Brightness

adjustment

: Manually drag to change the overall brightness of the screen.Often used in conjunction with contrast adjustment to achieve the best results.

Exposure         Exposure         Exposure	Automatic exposure	:	When checked, the device will automatically expose according to the environment used. When checked, this option is effective. When not checked, it can be adjusted by "exposure". Drag to adjust the exposure parameters of the device according to the use of the environment. This parameter is adjusted on
			the basis of automatic exposure. This parameter is only valid if the "automatic
			exposure" option is not checked.
Area Exposure	Regional exposure	:	Select an area as the reference for exposure.
Single Exposure	Single exposure	:	Perform once exposure.
R         -         B         -         -         -         -         -         -         -         -         -         -         -         -         -         -         -         -         -         -         -         -         -         -         -         -         -         -         -         -         -         -         -         -         -         -         -         -         -         -         -         -         -         -         -         -         -         -         -         -         -         -         -         -         - <td< td=""><td>Automatic white balance function</td><td>:</td><td>When checked, the device will perform automatic white balance according to the environment used;When not checked, it can be adjusted by "manual white balance" below.Note: abnormal color may occur when taking rich color images.A manual calibration of the automatic white balance is required.To do this, place a white paper on field of view, remove the "automatic white balance" box and check it again.Let the device perform an automatic white balance.When you're done,</td></td<>	Automatic white balance function	:	When checked, the device will perform automatic white balance according to the environment used;When not checked, it can be adjusted by "manual white balance" below.Note: abnormal color may occur when taking rich color images.A manual calibration of the automatic white balance is required.To do this, place a white paper on field of view, remove the "automatic white balance" box and check it again.Let the device perform an automatic white balance.When you're done,
K     +       B     -       -     +       112	Manual white balance function	:	unchecked the "auto white balance" box. Drag to adjust the white balance parameter of the device according to the usage environment. This parameter is adjusted up and down on the basis of automatic white

			balance. This parameter is only valid if the
			"automatic white balance" option is not
			checked.
Area White Balance     Imp	Areal white balance	:	Select an area as the reference for the white balance.
🖉 Single White Balance	Single white	:	Perform a white balance.
WB S	balance		
Contrast	Contrast adjustment	:	Manually drag the screen to change the
¥ ·	function		contrast between light and dark, often used in conjunction with brightness adjustment for
			best results.
Saturation	Saturation control	:	Manually drag to change the color saturation
	function		of the image for best results.
667 Sharpness 	Definition function	:	
Power Line	Power line	:	Note: the power frequency currently supports
50Hz 🗸	regulation		50/60hz, please choose according to the national power frequency.
	Restore default		Restore to the factory settings
Π	Restore default	•	
🖓 🗸 Save Settings	Save parameter	:	When checked, save the current settings when
			the device is power-off or the software is
			turned off.

# **3.4 Measuring tool**

m

1

1/

n

Calibration function and measurement function can accurately measure in editing area.

Note: calibration is the first step for accurate measurement.

	Ê	-0	×	
	Measurer	ment Tools		
	H-H	~ "//	/1/1	
		$\bigcirc \supset$	Ļ.	
	X,	<u> </u>	$\bigcirc$	
	$\bigcirc$	$\bigcirc$	⊕⊕	
	ଜକ ।			
	Calibratio			
	Unit mn	1	▼	
	Line type		<b>•</b>	
	Line widt	h		
		• +	8	
	Transpar — 🏾		0 %	
	Color			
Linear measurement	: M	leasure th	ne distar	nce as the crow flies between two points.
Parallel line measurement	: M	leasure th	ne distar	nce between parallel lines.
		r .1	1.	1
Double parallel line measurement	: M	leasure th	ie distai	nce between two parallel lines.
Stitche	: M	leasure th	ne circu	mference of the stitches.
Measurement				

Rectangular measurement : Measure the circumference and area of a rectangle.

$\bigcirc$	Polygon measurement	:	Measure the perimeter and area of a polygon.
$\mathcal{A}$	Arc measurement	:	Measure the length of the arc.
A	Four-angle measurement	:	Four points measure angles.
٢	Angle measurement	:	Measure angles at three points.
0	Oval measurement	:	Measure the radius, circumference, and area of an ellipse.
0	Circular measurement Circle measurement	:	Measure the area, circumference, and radius of a circle.
Θ	(circle center, radius to determine the circle)	:	Measure the area, circumference, and radius of a circle.
	Circle measurement (three points to determine the circle)	:	Measure the area, circumference, and radius of a circle.
Ō			
0	Circle measurement	:	Measure the circle tool.
₽₽	Circular measurement (concentric radius center distance)	:	Concentric radius center distance measuring tool.
ውብ	Circle measurement (three-point circle center distance)	:	Three-point circle center distance measuring tool.
ليتينا	Scale bar		
יויויו	Calibration	:	Set the calibration.Click this function to bring up the calibration Settings window, as shown in the figure

lame	Length	Pixel	Unit	Note	Name
efault	100	100	mm		
					Actual Length
					Pixel Length
					Unit
					mm IT
					Decimal Places
					3 💌
					Additional Remark
		III			•

calibration can be set in the window.

mm IV	unit	:	The measurement display unit of the current preview screen.
	Linetypes	:	Sets the drawing line shape. There are five linetypes to choose from.
+ 8	Line width	:	Sets the pixel value of the width of the drawing line, ranging from 1 to 10 pixels.
- [	Transparency	:	Sets the transparency of the drawn line.
	Color	:	Sets the color of the drawn line.

# **3.5 Drawing Tools**

According to the requirements, draw in the editing screen, as shown in the figure:

Drawing i	Drawing Tools 🛛 🗙						
	ø						
$\sim$	5	l					
_,	€_,	2					
≥,	$\diamond$						
		8					
$\triangle$	$\[\]$	$\bigcirc$					
Line Type							
Line Widtl — -{	י 	8					
Transparency - [] + 0 %							
Color							

	Pencil
<i>b</i>	Pen
$\mathbf{i}$	Straight Line
$\mathbf{M}$	Straight Line with Single Head Arrow

5	Straight Line with Double-headed	
2	Curve	
2	Arrow Curve	
1	Broken Line	
l,	Single Arrow Polyline	
, L	Double Arrow Polyline	
	Rectangle	
	Rounded Rectangle	
	Parallelogram	
$\diamond$	Rhombus	
8	Rhombus	
$\bigtriangleup$	Isosceles Triangle	
	Right Triangle	
	Cylinder	
Ø	Cube	
<b> </b> ▼	Line Type	Sets the drawing line shape. There are five linetypes to choose from.
+ 8	Line Width	Sets the pixel value of the width of the drawing line, ranging from 1 to 10 pixels.
- [	Transparency	Sets the opacity of the drawn line.
	Color	Sets the color of the drawn line.
3.6 Annotati	ion Tools	




Т	Input text.
A	Set font, size, effect.
	Sets the display color of the inserted font.
	Insert picture.
	Take a screenshot to the clipboard and draw the square area to the clipboard.
■ 截图时隐藏当前窗口	When ticked, the screenshot hides the current program, otherwise it displays
	normally.
٢	Take a screenshot of any drawn area.
	Take a screenshot of the rectangular area.
<u>.</u>	Take a screenshot of the entire window area.
34	Take a screenshot of the entire screen area.

### 3.7 Eraser

Erase the "object" in preview screen

As shown in the figure:



∎়	Rrase Objects	:	Erase the entire selected object.
	Rrase All	:	Erase all objects on screen.

# **3.8 Image Filters**

**k**ohd

Monochrome

 $\star$ Image Filters only available for static images

Operate the image in the editing area according to the requirements, as shown in the figure:



₩ <u></u>	Proportional Crop	:	Cut images in equal proportions Click and the window on the right will appear	
Ę	Resize image	:	Adjust the size of the image. Click and the window on the right will appear	Resize Image         Resize Image         Basis ( B )       Percentage       Resolution         Horizontal (H)       2592         Vertical (V)       1944         Maintain Aspect Ratio (M)         Confirm       Cancel
	Inverted color			

ion





### 3.9 Calibration and Measurement

Step 1: Acquire benchmark image.

Place a scale ruler in the filed of view. Please use the same conditions (i.e. same objective lens, same camera resolution) as the target observation ojects. If a magnification reduction lens is used to measure the target observation ojects, then a magnification reduction is also needed to measure the scale bar. This is what we called acquire benchmark image.

As figure follows:




Step 2: Click "Calibration" in the left measuremtn tools toolbar, calibration table will shown as below figure.

Calibration Tai	ble				- ×
Name	Length	Pixel	Unit	Note	Name
default	100	100	mm		
					Actual Length
					Pixel Length
					Unit
					Decimal Places
					3 💌
					Additional Remark
•		111			•
	Apply	Clo	se		Add Edit Delete

Step 3: Setting scale ruler.

Draw a straight line on the scale ruler of benchmark image, pixel value will be automatically updated at "pixel length" on calibration table. Shown as below figure :



ame	Length	Pixel	Unit	Note	Name
х	10	1664	μm		10X
					Actual Length
					10
					Pixel Length
					1664
		<u> </u>			Unit
					μm 💌
					Decimal Places
					3 🔻
					Additional Remark
			_		
					>

When drawing the straight line, the longer scale length, the more accurate results. For example, 10 scale lengths accuracy will be better than 1 scale length.

Step 4: Input the actual length and unit of the ruler in the calibration table. As shown in the figure below:

L			10X	
Actu	al Leng	yth		
	10	)		
Pixel	Lengt	n		
	16	54		
Unit				
	μm			
Decir	mal Pla	ices		
	3			
Addit	ional F	Rema		
Г				

Step 5: click "Add" to save to the calibration table.

Calibration To	able				- ×
Name	Length	Pixel	Unit	Note	Name
10X	10	1664	μm		10X
					Actual Length
					10
					Pixel Length 1664
					Unit
					μm  ▼
					Decimal Places
					3 🔍
					Additional Remark
<					
	Apply	Clo	ose		Add Edit Delete

lame	Length	Pixel	Unit	Note	Name		
0X	10	1664	μm		- Nume	10X	
					Actual Leng		
i i					10		
					Pixel Length		
					166	4	
		-			Unit	-	
	 	_	2		Decimal Pla	<b>•</b>	
					3	l▼	
					Additional R		
					> 1	2	3


1)"Add": Add the calibraion information to the list on the left

2)"Edit": Revise or edit selected calibration item on the left

3)"Delete": Delete selected calibration item on the left

4) "Apply": Activate selected calibration item as benchmark for current image


Step 6: After activating selected calibration item , put the object into the field of view, and click " " in the left bar.Select the corresponding measurement tool, and the measurement results will be displayed in the right display area. As shown in the figure below:




1.Polygon Measurement Perimeter=1478.643 mn Area=118422.500 mm <sup>2</sup>	n
2.Straight Line Distance Measurement Length=555.306 mm	
3.Parallel Lines Distance Measurement Length=284.580 mm	

Step 7: Export the measurement results.

Click in the display bar on the right to export the data to Word or Excel.Enter the file name and path, and click "save". As shown in the figures below:





## 3.10 Thumbnail display area

General status



Right click the mouse on the thumbnail




- Open:Open the image in the system default mode
- Open Directory:Open the directory where the selected image is located
- Copy:Copy the currently selected image
- Paste:Paste a copied image (note: only paste in thumbnails)

- Delete:Delete the selected image and delete the source file
- Rename:Rename the selected image
- Send to Email:Add the selected image as an attachment to OUTLOOK
- Refresh: Refresh the thumbnail area of the current image

Click the bottom button to switch to display the

measurement data



All	▼
-----	---

Picture display type

All, Photo, PDF and Video are options.

All	▼
All	_
Photo	
PDF	
Video	

Thumbnail

:

Display pictures and files under the specified path.

17	Measurement	:	Display measurement data	
	data			
W	Export to Word	:	Export the currently displayed	
			measurements to Microsoft Office	
			Word	
XI	Export to	:	Export the currently displayed mea	asurements to
7	EXCEL		Microsoft Office Excel	