# Brugsvejledning for 2648.00 benzinkanon

25.04.12

AC 2648.00

## Benzinkanon

Kanonen består af følgende dele:

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AC 2648.00

### Gasoline cannon

The cannon consists of the following parts:

- A strong transparent plastic tube with an opening in both ends.
- A piezoelectric igniter (97700630).
- A rubber stopper with hole (0455.10).
- A cork (97570150).

The igniter shall be mounted in the rubber stopper with the large diameter end towards the igniter handle. To reduce friction, a small amount of liquid soap can be applied.

Place the stopper with igniter in one end of the tube (see Fig.1).

Fig. 1

### Firing of the cannon

A cotton cloth or similar wetted with benzine (petroleum spirits, ligroin) and brought into the plastic tube (see Fig. 2). Let the cloth stay in the tube for 10 to 15 seconds then drag out the cloth and insert the cork loosely in the open end of the plastic tube.

Fig. 2

The cannon can now be fired. It may require a few sparks from the igniter before it goes off.

### Warning!

It is recommended to secure the cannon with one or two experimental stands. Ensure that spectators are not in the line of fire.

The operator (and students) should wear protective glasses. It is recommended that a Perspex protection screen is placed between the cannon and the spectators.

The cork must be pressed only loosely into the pipe to avoid the risk of bursting the pipe.

Use benzine for the experiment.

Do not under any circumstances use alcohol in the tube since this can cause cracks and other defects in the acrylic tube which may fracture when the cannon is fired.

### Note!

Before using the equipment, check the tube for cracks. Even very small cracks can cause the tube to fracture when the cannon is fired.

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© A/S Søren Frederiksen, Ølgod