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Guide for pH electrodes

1. Storage of the electrode:

The pH membrane (the lower round part of the electrode) should always be kept moist. For storage, the protective cap must contain 3 M potassium chloride or possibly buffer pH 3-4. Storage bottle 541600 can be used with advantage.

2. The pH-sensitive membrane should not be touched

Damage such as scratching or scraping will result in a slow and inaccurate response.

3. Check that the internal liquid completely fills the pH-sensitive membrane

Release any air bubbles by shaking the electrode in a downward direction.

4. For normal use, rinse the electrode with demineralised water. Clean off grease with acetone

5. The pH meter can be checked as follows:

Short the pH connector with clips or shorted BNC connector. With the buffer button in the centre position, the pH reading should be approx. 7. Furthermore, the range should be adjustable from approx. 5.7 to 8.3 pH.

6. Electrode lifespan:

If the buffer knob or "slop screw" does not adjust the pH reading to the correct value, check the pH meter. If this works satisfactorily, clean and soak the electrode in 0.1 N HCl for 2 hours, then place it in buffer pH 7 for 1 hour.

If the adjustment is still unsuccessful and the buffer is in order, the electrode must be replaced.