2024-08-22 / DTS

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Anvendelse

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Application

Cs-137/Ba-137-Isotope Generator (T), (R) from Spectrum Techniques for extracting radioactive preparation for experiments with measurements of half-life. Approx. 1000 experimental doses of Ba-137 can be generated. The decay is a gamma decay.

When a small amount of Ba-137 has been milked from the generator, there is a short period of time when there is no significant amount of Ba-137 left in the isotope generator. You must therefore wait approx.. 15-20 minutes before a new equilibrium has been established in the generator.

The generator can be milked for isotopes 2-3 times in an hour, after which it must have time to 'recover'.

Materials

A case for storage.

The isotope generator with radioactive sources.

Disposable syringe used to supply the generator with extraction fluid.

10 small metal trays to collect the fluid that is milked from the generator and contains the radioactive isotopes.

Bottle with 200 mL of extraction fluid.

How to use

- 2 ml. of extraction fluid is drawn into the syringe.
- One of the small metal trays is put on the table ready to receive the radioactive isotope fluid.
- The protective caps are removed from the inlet and outlet of the isotope generator.
- The syringe is attached to the inlet (the long protrusion located at the top of the generator when the writing on the generator is facing correctly).
- The extraction solution is slowly injected into the generator. Remember to place the outlet over the small metal tray so the 'milked' fluid can be collected.
- Stop injecting extraction fluid into the generator when 7 drops of radioactive 'milk' have been collected on the metal tray.
- 7 drops of isotope 'milk' will measure as an activity of 300 Becquerel (Bq) just above the metal tray.
- The measurement of radioactive decay can now begin.

- After 15 min, the activity is approx. 4 Bq and must be discarded into the sink. (The half-life of Ba-137 is 2.6 min.). Remember to rinse afterwards.

Security measures

Make sure that the syringe is correctly mounted on the isotope generator – i.e. the inlet port.

Make sure not to draw the liquid back into the syringe, as this can damage the non-return valve in the generator.

Only use original extraction fluid.


Store the isotope generator in a locked cabinet.

Disposal of the generator must be done according to laws and regulations.

The warranty does not cover damage to the isotope generator caused by misuse or the use of non-original extraction fluids.

Theory

The isotope decay is illustrated in the diagram below:

Restitution time for the isotope generator:

Minutter	Procent
1,0	Ca. 23,0%
2,6	Ca. 50,0%
5,2	Ca. 75,0%
12,0	Ca. 96,0%
15,0	Ca. 98,0%
20,0	Ca. 99,9%

The half-life of Ba-137* is 2.6 minutes, and the energy of the gamma quanta is 662 keV.