## Radioaktive sten 670255

English  
manual  
on p. 3

2016.06.02 / HS

AC 670255

### Beskrivelse

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AC 670255

### Description

This set consists of samples of radioactive minerals and must be stored out of reach of children.

Note that the rocks may be brittle after hundreds of millions of years of exposure to radiation – handle them carefully. Please wash your hands after working with these samples.

The rocks can be utilised in many different experimental contexts, from a simple phenomenological demonstration of radioactivity to an actual identification of some of the constituent isotopes.

The elements has their natural distribution of isotopes and as the minerals are very old, all daughter nuclei in the decay chains will be present in equilibrium.

### Contents

The mineral samples are labelled with numbers. (A more detailed description can be found overleaf.)

- 1 - Uraninite
- 2 - Thorianite
- 3 - Euxenite
- 4 - Allanite
- 5 - Pyrochlore
- 6 - Rødberg

### Suitable equipment

512515	GM tube on post, BNC plug
512525	GM tube large area, BNC plug
513600	Geiger counter
513570	GM sensor, jack plug (For e.g. Pasco dataloggers)
518500	Scintillation detector for 518000
518000	Multichannel analyser

## **1 - Uraninite**

Uraninite (pitchblende),  $\text{UO}_2$ , is the most uranium rich mineral with up to 88 % uranium. It also contains between 5 and 20 % lead and helium of radioactive origin, together with radium and varying traces of rare earths. Uraninite almost always occurs associated with granitic rocks. In Norway, uraninite often occurs in so-called granite pegmatites. Most pegmatites are associated with granites and consist like these of coarse-grained feldspar, quartz and mica. Here, natural enrichment processes has concentrated more uranium, thorium and other parts than in ordinary granite. By running a Geiger-Muller tube across the specimen, one will find one or more greenish black areas where you get higher counts. This is uraninite which often has thin layers of greyish yellow, later formed uranium minerals close by. On a larger scale the same method can be used to find uraninite or other radioactive minerals in nature.

## **2 - Thorite**

Thorite,  $\text{ThSiO}_4$ , is a thorium rich mineral with up to 82 % thorium. Thorite occurs like uraninite partly in granite pegmatites and partly in corresponding coarse formations as nefelin syenite pegmatites from quarts-poor rocks. In the sample, thorite can be seen as a glossy black to matte brown, very porous mineral.

## **3 - Euxenite**

Euxenite ( $\text{Y, Ca, Ce, U, Th})(\text{Nb, Ta, Ti})_2\text{O}_6$ ) is a glossy black mineral that often occurs in Norwegian granite pegmatites. Euxenite may occur in lumps of up to several kilograms. They usually contain significant amounts of uranium and thorium (up to 10 %) and will in larger pieces provide substantial count rates with a Geiger counter. The mineral is also of interest in terms of yttrium, niobium and tantalum, which have an ever increasing use in electronic components.

## **4 - Allanite**

Allanite (Orthit) ( $\text{Ce, Ca, Y}_2(\text{Al, Fe}^{3+})_3(\text{OH})$ ) is one of the most common minerals in the Norwegian pegmatite deposits. They can occur in very large crystals of over 100 kg. They have a variable composition including 0 to 3% uranium and thorium. The mineral can also contain yttrium and cerium.

## **5 - Pyrochlore**

Pyrochlore ( $\text{Na,Ca}_2\text{Nb}_2\text{O}_6(\text{OH},\text{F})$ ) can contain a few percents of yttrium, cerium, thorium and uranium. The mineral can be seen as small brown to black vitreous grains. Sometimes it is seen as larger deposits and it is then interesting in relation to extraction of niobium and the others elements mentioned. This sample is from nefelin syenite pegmatites in the Langesund area where it occurs as small grains along with other rare minerals.

## **6 - Rødberg**

Rødberg (raudberg) consists primarily of dark red siderite,  $\text{FeCO}_3$ , together with some metallic hematites,  $\text{Fe}_2\text{O}_3$ . This type occurs at Fen, Ulefoss in Telemark where it, together with other carbonate-containing rocks is formed by solidification of melts in the crust of the earth. This rare and special mineral was previously prospected due to the content of thorium and rare earths. The sample is only weakly radioactive and sufficiently long counting times must be used, preferably with an extra sensitive GM tube.

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© Frederiksen Scientific A/S