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hentes på vores hjemmeside.

Air Track

The air track is a precision instrument for physics education, where an airflow reduces friction between the track and the carts. This makes it possible to perform experiments with high accuracy and obtain stable, reproducible measurement results.

The track is made of aluminium with an anodized surface, which protects against wear and ensures a long service life. A millimetre scale is mounted on one side of the track for precise distance measurements.

Scope of Delivery

The air track is supplied with:

- 1 air track, 150 cm
- 2 carts
- Standard accessory set in 2 plastic boxes
- Wide and small support foot

Technical Specifications

- Length: 150 cm
- Width: 85 mm (widest point)
- Height: 115 mm (base to top)
- Material: Anodized aluminium
- Height adjustment: 3 adjustable feet
- Scale: Millimetre scale on one side

Operation

Connect the hose to the air inlet of the air track, then connect the blower to a power outlet using the mains cable. The airflow is adjusted using the stepless control knob.

Optimal airflow is achieved when the carts are just suspended above the track without touching it. Insufficient airflow may cause friction between the carts and the track, while excessive airflow may create turbulence around the carts.

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© Frederiksen Scientific A/S, Ølgod.