

COMPATIBILITY CHART

Description	Cat. No.	CAPP			Gilson			Eppendorf			Brand			Thermo			Sartorius			HTL		
		Bravo / scopiPipette / CAPP Aero / CAPP Maestro / CAPP Tronic	single and multichannel pipettes	0.2-2 µL 2-20 µL 5-50 µL 10-100 µL 20-200 µL 100-1000 µL 1000-5000 µL 1000-10000 µL	Pipetman / Pipetman G	single and multichannel pipettes	0.2-2 µL 0.5-10 µL 2-20 µL 10-100 µL 20-200 µL 100-1000 µL 1000-10000 µL	Research Plus	single and multichannel pipettes	0.2-2 µL 0.5-10 µL 2-20 µL (grey) 2-20 µL (yellow) 10-100 µL 20-200 µL 100-1000 µL	TransferPipette, TransferPipette S	single and multichannel pipettes	0.5-10 µL 2-20 µL 10-100 µL 20-200 µL 100-1000 µL	Fimmpipette FL / Novus	manual single and multichannel and electronic single and multichannel pipettes	0.5-10 µL 2-20 µL 10-100 µL 20-200 µL 100-1000 µL 1000-10000 µL	Proline Plus / Picus	manual single and multichannel pipettes and electronic single and multichannel pipettes	0.2-10 µL 0.5-10 µL 5-120 µL 10-300 µL 50-1000 µL	Labmate Pro (same as Labnet BioPipette)	manual single and multichannel pipettes	0.2-2 µL 2-20 µL 10-100 µL 20-200 µL 100-1000 µL 1000-5000 µL 1000-10000 µL
CAPP Expell Filter Tips																						
Expell 10 µL, w/ filter, clear, sterile	5130030C																					
Expell 10 µL, w/ filter, clear	5131030C																					
Expell 10 µL XL, w/ filter, clear, sterile	5130060C																					
Expell 10 µL XL, w/ filter, clear	5131060C																					
Expell 10/20 µL XL, extra narrow, w/ filter, clear, sterile	5130061C																					
Expell 10/20 µL XL, extra narrow, w/ filter, clear	5131061C																					
Expell 20 µL, w/ filter, clear, sterile	5130062C																					
Expell 20 µL, w/ filter, clear	5131062C																					
Expell 100 µL, w/ filter, clear, sterile	5130066C																					
Expell 100 µL, w/ filter, clear	5131066C																					
Expell 200 µL, w/ filter, clear, sterile	5130090C																					
Expell 200 µL, w/ filter, clear	5131090C																					
Expell 1000 µL, w/ filter, clear, sterile	5130150C																					
Expell 1000 µL, w/ filter, clear	5131150C																					
Expell 1250 µL XL, w/ filter, clear, sterile	5130123C																					
CAPP Expell and ExpellPlus individually wrapped Pipette Tips																						
ExpellPlus 10 µL, clear, sterile, indiv. wrapped	5030012C																					
ExpellPlus 10 µL XL, clear, sterile, indiv. wrapped	5030042C																					
ExpellPlus 200 µL, clear, sterile, indiv. wrapped	5030072C																					
Expell 1000 µL (1250 µL), clear, sterile, indiv. wrapped	5130132C																					
CAPP Expell and ExpellPlus Filter Tips in PaperBox Racks																						
ExpellPlus 10 µL, w/ filter, clear, sterile	5030029C																					
ExpellPlus 10/20 µL XL, extra narrow, w/ filter, clear, sterile	5030059C																					
ExpellPlus 20 µL, w/ filter, clear, sterile	5030063C																					
ExpellPlus 100 µL, w/ filter, clear, sterile	5030064C																					
ExpellPlus 200 µL, w/ filter, clear, sterile	5030089C																					
ExpellPlus 300 µL, w/ filter, clear, sterile	5030114C																					
Expell 1000 µL, w/ filter, clear, sterile	5130149C																					
ExpellPlus 1250 µL XL, w/ filter, clear, sterile	5130119C																					
CAPP Expell Pipette Tips Reload System																						
Expell 10 µL, clear, reload	5030035C																					
Expell 10 µL XL, extra narrow, clear, reload	5030055C																					
Expell 200 µL, clear, reload	5030095C																					
Expell 300 µL, clear, reload	5130111C																					
Expell 1000 µL (1250 µL), clear, reload	5130155C																					

This compatibility chart is for your reference. We are not responsible for technical modifications or changes made to the mentioned pipettes other than those of CAPP. Specifications like inaccuracy and imprecision values have not been checked. For a complete compatibility list please contact info@cappahn.com.

■ compatible x not compatible