

Elektrodeholder, dobbelt, universal, 449710



Anvendelse

Elektrodeholder til to elektroder. Kan benyttes til både runde elektroder og pladeformede elektroder. Elektrodeholderen er designet således den kan hvile på kanten af et 250 mL bægerglas, lav form.

Tilslutning

Tilslutning udføres med 4 mm sikkerhedsstik.



Elektrodetyper

- *Pladeelektrode*
 - 449800 - Pladeelektrode, kobber, 50 x 87 mm
 - 449805 - Pladeelektrode, aluminium, 50 x 87 mm
 - 449810 - Pladeelektrode, zink, 50 x 87 mm
 - 449830 - Pladeelektrode, jern, 50 x 87 mm
- *Stang elektroder (150 mm)*
 - 451000 - Kobberelektrode, 6 mm
 - 451010 - Nikkelelektrode, 6 mm
 - 451020 - Aluminiumelektrode, 6 mm
 - 451030 - Tinelektrode, 6 mm
 - 451040 - Jernelektrode, 6 mm
 - 451050 - Kulelektrode, (250 mm), Ø6 mm
 - 451060 - Zinkelektrode, 6 mm

Universal electrode holder, double, 449710



Operation

Electrode holder for two electrodes. Can be used for both round electrodes and plate-shaped electrodes. The electrode holder is designed so that it can rest on the edge of a 250 mL beaker, low form.

Connections

Connection is made with a 4 mm safety plug.



Electrode types

- *Plate electrode*
 - 449800 - Plate electrode, copper, 50 x 87mm
 - 449805 - Plate electrode, aluminium, 50 x 87mm
 - 449810 - Plate electrode, zinc, 50 x 87mm
 - 449830 - Plate electrode, iron, 50 x 87mm
- *Round electrode (150mm)*
 - 451000 - Round electrode, copper, 6mm
 - 451010 - Round electrode, nickel, 6mm
 - 451020 - Round electrode, aluminium, 6mm
 - 451030 - Round electrode, tin, 6mm
 - 451040 - Round electrode, iron, 6mm
 - 451050 - Round electrode, (250mm) carbon, 6mm
 - 451060 - Round electrode, zinc, 6mm